

UNIT - Ist

UV VISIBLE SPECTROSCOPY

- Analysis →

It is the quantitative & qualitative determination of any sample.

- Spectroscopy →

$\text{Spectrum} + \text{scopy}$
 light
 to determine


• It is the branch of science which is used for quantitative analysis of matter / sample by using light (Electro-magnetic radiations) --

eg. amount of substance in sample etc...

- $\boxed{\text{EMR}} \rightarrow \text{Electromagnetic Radiations}$

- light flows in the form of wave

- It contains electric & magnetic fields perpendicular to each other --



- Wavelength \rightarrow Distance b/w the two consecutive waves (crest/trough). $\rightarrow \lambda$

- frequency \rightarrow No. of waves (cycles) in one sec. \rightarrow hertz (Hz) ..

frequency, $\nu = \frac{c}{\lambda}$ - velocity of light
 λ - lambda

• Gamma Rays	< 0.001 nm
• X-rays	< 0.01 - 10 nm
• U.V	200 - 400 nm
• Visible light	400 - 800 nm
• Infra Red	0.8 - 200 μ m
• Microwaves	0.01 - 1m
• Radio waves	1 - 10^7 m

↑
E N E R G Y ↓

• UV Visible Spectroscopy →

Used for quantitative analysis of matter/sample by using UV (ultra violet) and Visible light.

• the overall range of wavelength of UV + visible is 200 - 800 nm.

• Principle →

- EMR contain energy (photons)

$$E = h\nu$$

∴ $h = \text{Planck's constant}$

$$6.626 \times 10^{-27} \text{ erg} \cdot \text{sec}$$

$$\nu = \frac{c}{\lambda}$$

$$E = h \cdot \frac{c}{\lambda}$$

$$E \propto \frac{1}{\lambda}$$


Frequency ↑ = Wavelength ↓

Wavelength ↓ = Energy ↑



When light (EMR) is passes through a sample

Principle →


- EMR light contain packets of energy (photons) -
- Now, when this light passes through the sample, \rightarrow particles/molecules of sample absorb this light.



So, particles get excited -




- Now, UV-Visible spectroscopy's principle depends upon the absorption of light.



because particles absorb some light.

eg. — who absorb more light — ??



- So, if we determine the intensity of light, then we calculate the concentration of particles —

Absorption \propto concentration of particles

Absorption \uparrow = concentration \uparrow


So, Now, we have to calculate the

Absorption

How — ??


eg.

Alkali' delay



principle based on

Beer-Lambert law



Acc to this,


light ~~A~~ Absorption is directly proportional to path length and the concentration of the solution.

$$A = \epsilon \times l \times c$$

ELECTRONIC TRANSITIONS

The absorption of electromagnetic radiation causes electrons to be excited, which results in promotion from a bonding / non-bonding orbitals to an anti-bonding orbitals i.e. $n \rightarrow \pi^*$, $\pi \rightarrow \pi^*$

OR



• three types of electrons involving in organic molecules :-

i) σ -electrons :- these electrons are involved in σ -bonds, these

are the lowest-energy occupying molecular orbital i.e. σ -orbitals / σ -bonds

ii) π -electrons \rightarrow these electrons are involved in π -bonds in unsaturated hydrocarbons.


these π -orbitals / π -bonds lie at higher energy levels.

iii) n -electrons \rightarrow these electrons are involved in non-bonding orbitals in lone pair or non-bonding (n) orbitals

- they also lie at higher energy level

(eg) Nitrogen, Oxygen & halogens etc.

The possible electronic transitions can graphically shown as :-



• The possible electronic transitions are :-

- i) $\sigma \rightarrow \sigma^*$ transition
 - ii) $\pi \rightarrow \pi^*$ transition
 - iii) $n \rightarrow \sigma^*$ transition
 - iv) $n \rightarrow \pi^*$ transition
 - v) $\sigma \rightarrow \pi^*$ transition
 - vi) $\pi \rightarrow \sigma^*$ transition
- } theoretical possible
(show absorption in region above 200 nm)

• The most probable transition is from the highest occupied molecular orbital (HOMO) to the lowest unoccupied molecular orbital (LUMO).

- i) $\sigma \rightarrow \sigma^*$ transition :- A transition of an electron from a sigma bonding orbital to the higher energy anti-bonding sigma orbital ($\sigma \rightarrow \sigma^*$)

• high energy required $\left[\because \frac{1}{\lambda} = E \right]$

• very short wavelength



show wavelength of 125 nm

- ii) $\pi \rightarrow \pi^*$ transition \rightarrow π electron in a bonding orbital is excited to

corresponding anti-bonding π^* orbital.

• Available in compounds with unsaturated centres.

(eg) Alkenes, —HC=CH—CH=CH—
alkynes, carbonyl, etc.

Alkenes generally absorb in the region
170 to 205 nm.

iii) $n \rightarrow \sigma^*$ transition \rightarrow It involves saturated

compounds with atoms containing

low pair of electrons O, N, S

(eg) NH_3



wavelength \rightarrow 150 - 250 nm.

iv) $n \rightarrow \pi^*$ transition \rightarrow It involves an electron


of atoms containing low pair is
excited to π^* anti-bonding orbitals

Required least energy & longer wavelength


(eg) C=O: , $\text{HC}\equiv\text{N:}$ etc.

wavelength \rightarrow 280 - 300 nm

Energy :-



(eg)



Three types of
electrons


CHROMOPHORES

Chromophores are the molecules or parts of molecules that absorb light in the UV-Visible region and responsible for imparting colour to the compounds.

→ reflect a certain color.

(eg) Nitro-compound ($-NO$) are generally yellow in color.

Colorless benzene



→ Pale Yellow color

+ hydroxyl group → dark yellow.

Autochrome

• Some of the important chromophores are

ethylenic, acetylenic, Carbonyl, ester,

Nitrils group etc.

$\pi-\pi$, $n-\pi$

Types

Independent chromophores

• A single chromophore imparts color to the compound.

Dependent chromophores

• More than one chromophore is required to produce color in compound.

(eg) azo group ($-N=N-$)

Nitro group ($-NO$)

(eg) $C=O$ groups,

$C=C$ groups etc.

Acetone

one ketonic group → colorless

2 ketonic group → Yellow


3 ketonic group → Orange

Auxochrome

It is a group of atoms that get attached to the chromophore and increase the colorfulness of the chromophore.

• Also called as Color enhancing groups.


eg.



Benzene

 $\lambda_{\text{max}} \rightarrow 255 \text{ nm}$ $\epsilon_{\text{max}} \rightarrow 203$

Absorption maximum



Aniline

 $\lambda_{\text{max}} \rightarrow 280 \text{ nm}$ $\epsilon_{\text{max}} \rightarrow 1430 \text{ nm}$

this group act as a chromophore.

- $-\text{OH}$, $-\text{OR}$, $-\text{NH}_2$, $-\text{NHR}$, $-\text{NR}_2$, $-\text{SH}$ etc.

• Also auxochromic groups

increases the wavelength of molecules.

Mechanism \rightarrow


they contain non-bonding electrons, which

extended the conjugation of the

chromophores by sharing the non-bonding electron.

SPECTRAL SHIFTS

The position of absorption maximum and absorption intensity can be modified in different ways by some structural changes or solvent changes :-



i) Bathochromic shift \rightarrow It is an effect in which (Red shift) the absorption maximum (λ_{max})

is shifted towards longer wavelength due to presence of an auxochrome or by the

change of solvent.

eg. The $n \rightarrow \pi^*$ transition for carbonyl compounds experiences bathochromic shift (red shifts) when the solvent polarity is decreased.

ii) Hypsochromic shift \rightarrow (Blue shift)

It is an effect in which the absorption maximum is shifted towards shorter wavelength.

iii) Hypsochromic shift \rightarrow

It is an effect in which

the intensity of absorption maximum increases i.e. the value of ϵ_{max}

eg. Auxochrome increases the intensity of Absorption.

iv) Hypochromic shift \rightarrow

It is an effect in which the intensity of absorption maximum decreases i.e.

the value of ϵ_{max} decreases.

It may be caused by the introduction of group which distorts the geometry of molecules.

eg. Biphenyl - ϵ_{max} 19000 whereas 2-methyl biphenyl at ϵ_{max} 10250

SOLVENT EFFECT ON ABSORPTION SPECTRA

The choice of solvent used in UV spectroscopy is most important because solvent can change the absorption max.

• A most suitable solvent is one that does not itself absorb radiation in the region under investigation (UV-VIS radiation).

Polarity \downarrow \rightarrow Absorption maximum (λ_{max})

wavelength \uparrow

Solvent polarity

- Most commonly used solvent is 95% ethanol.

Solvents	λ of Absorption
Ethanol	210 nm
Hexane	210
Methanol	210
Cyclohexane	210
Water	205
Benzene	280
Chloroform	245

← Polarity

Water > Methanol > Ethanol > Benzene > Hexane

- The position as well as the intensity of absorption maximum get shifted by changing the solvent polarity.

- eg. The α, β -unsaturated carbonyl compounds shows two different shifts.

i) $n \rightarrow \pi^*$ → Absorption band moves to shorter wavelength by increasing the solvent polarity

• λ_{max} of acetone is 279 nm in hexane as compared to 264 nm in water

ii) $\pi \rightarrow \pi^*$ → Absorption band moves to longer wavelength by increasing the solvent polarity

iii) $n \rightarrow \sigma^*$ → Same as $n \rightarrow \pi^*$ (shorter wavelength)

BEER & LAMBERT'S LAW

- Beer's law \rightarrow Absorbance is directly proportional to the concⁿ of solⁿ

$$A \propto c$$

- Lambert's law \rightarrow Absorbance is directly proportional to the length of solⁿ

$$A \propto l$$

On combining,

$$A \propto cl$$

$$A = \epsilon cl$$

where, ϵ = molar absorption coefficient


c = concⁿ of solⁿ (mol/litre)

l = length of light pass through sample

* Derivation \rightarrow

Acc to Beer-Lambert's law

When a monochromatic light is passed through a solution containing the absorbing substance, the decrease in the intensity of light with path length is proportional to the concentration of solution and the intensity of light.



$$-\frac{dI}{d\lambda} = KIC$$

where,

I = intensity of light

c = concn (mol/litre)

l = thickness

K = constant

Interchange,

$$-\frac{dI}{I} = K \cdot C \cdot d\lambda$$

Integrating both sides,

$$-\int_{I=I_0}^{I=I} \frac{dI}{I} = KC \int_{\lambda=0}^{\lambda=l} d\lambda$$

$$-\ln \frac{I}{I_0} = KC l \left[\ln \frac{I_0}{I} = KC l \right]$$

$\ln(x) = 2.303 \log(x)$

$$2.303 \times \ln \frac{I_0}{I} = 2.303 \times KC l$$

$$\log_{10} \frac{I_0}{I} = \epsilon c l \quad \text{--- (1)}$$

$$T = \frac{I}{I_0}$$

Transmittance \rightarrow amount of light emitted after absorption.



$$A > 0 \\ T < 100\%$$

Now,

$$A = \log_{10} \left(\frac{I_0}{I} \right) \quad \text{--- (2)}$$

Rough condition $A \propto \frac{1}{T}$

put the value of T into eq. 2

$$A = \log_{10} \left(\frac{I_0}{I} \right)$$

$$A = \log_{10} \left(\frac{I_0}{I} \right) \quad \text{--- (3)}$$

acc to eq. (1) + (3)

$$A = \epsilon \cdot c \cdot l$$

\rightarrow Molar absorptivity (molar extinction coefficient)


Absorption

concentration + path length

Deviation →

A/c to law,

Absorption is directly proportional to concn of solution.



Types :-

- i) Real deviation / True deviation
- ii) Chemical deviation
- iii) Instrumental deviation / Spectral deviation.

i) Real deviation → It is due to higher concentration of molecules in solution.

Also due to refractive index of absorbing medium.

ii) Chemical deviation →

It is due to chemical


changes such as association, dissociation, pH changes etc. in absorbing medium.

(eg.)

phenol → 270 nm, phenoxide anion → 287 nm

iii) Instrumental deviation →

It is due to polychromatic radiation, it leads to negative deviations. Monochromators are used to prevent this.



INSTRUMENTATION

in UV-Visible Spectroscopy

- Source of radiation, wavelength selector, sample cells, detectors - photo tube, photomultiplier tube, photo voltaic cells, silicon photodiode.


Spectrophotometer :-

→ spectrophotometer is made up of two instruments :-

Spectrophotometer

Spectrometer photometer

produce light of any wavelength
used to measure the intensity of light.



The essential parts of a spectrophotometer are :-

- i) Source of Radiation
- ii) wavelength selector (Monochromator)
- iii) Sample cells or Cuvettes
- iv) Detector
- v) Recording system

i) Source of Radiation →

The best source light is the one which is more stable, more intense and which gives the range of spectrum

from 180 - 700 nm :

⊙ Hydrogen discharge lamp → In these lamps,

hydrogen gas is stored under

high pressure. When an electric discharge is passed through the lamp, excited hydrogen molecules

will be produced which emit UV radiation.

⊙ Deuterium lamp → It is similar to hydrogen discharge lamp, but filled with

deuterium (D_2) in place of hydrogen.

It offers 3-5 times more intensity than hydrogen lamp.

⊙ Xenon discharge lamp →

In this lamp, xenon at 10-30

atm pressure is filled ~~with~~ it and it has two tungsten electrodes.

- Greater intensity than hydrogen lamp.

⊙ Mercury Arc →

This contains mercury vapours but not continuous so not widely used.


ii) wavelength selector →

and slits.

It consists of monochromator

⊙ Monochromators → It is used to disperse the radiation

according to wavelength.





- The essential elements of a monochromator are an entrance slit, a dispersing element, and an exit slit.

- The dispersing element may be a prism or grating.


generally made of glass, quartz or fused silica.
widely used in UV spectrophotometers.

③ Prism Monochromator → two types →



③ Grating monochromator →

It provides an alternative means of producing monochromatic light.



It consists of series of parallel lines (grooves) which reflect through highly polished surface of glass, quartz or alloy halides (grating)...

o Slit →


Entrance slit for entry of light and exit slit passed the single radiation of desired wavelength.

* Wavelength selector are the device which get the light from source of radiation and passes the light with desired wavelength which is required for detection to the sample (cuvette).

(iii) Sample cells / cuvettes →

These are sample containers, which holds the liquid samples.

- they are transparent & usually made up with quartz or fused silica also silicate glass..



iv) Detectors →

Detectors used in UV-Visible spectrophotometer can be called as photometric detectors.

- These are those devices which converts light source into electrical signal which finally displayed by recording system.

• most commonly used detectors are:-

- photo tubes / photo-emissive cells
- photomultiplier Tubes (PMT)
- photovoltaic cells or barrier-layer cells


o Phototubes or photo-emissive cells →

- It is composed of an evacuated glass tube, which consist of a photocathode and a collector anode.


- The photocathode is coated with element of

high atomic volume like caesium, potassium or silver oxide, which can liberate (emit) electron when light radiation falls on it.

- This flow of electrons toward anode produce a current proportional to the intensity of light radiation.



c) Photomultiplier tubes (PMT) →



- The principle employed in this detector is multiplication of photoelectrons by secondary emission of electrons.
 - It consists of a photocathode and a series of anode (dynodes).
 - Each dynode is maintained at 75-100V higher than the preceding one & electron emission is multiplied by a factor
- Phototubes have better sensitivity when compared to photovoltaic cells and hence are widely used.

of 4 to 5 due to secondary emission of electrons, hence an overall factor of 10^5 is achieved.

PMT can detect very weak signals.

Now, finally the large number of electrons arrive at collector (anode).

The number of electron falling on the collector measures the intensity of light incident on the cathode surface.

The response time of PMT is 10^{-9} s.


9) Photovoltaic cells (barrier-layer cells)

- Also known as photoelectric cell and operates without the use of a battery

- It consists of a metal base plate (iron/aluminium) that act as an electrode, on this thin layer of semiconductor metal

(like selenium) is deposited.

- then this covered by a very thin layer of silver or gold that act as second collector electrode.



- When light-radiation falls on selenium layer, electron are generated at the selenium-silver interface.

- Thus electrons are collected by the silver, this creates a electric voltage difference b/w the two electrode (silver & base cell) and causes the flow of current.


- Now, this flow is directly proportional to the intensity of the incident radiation beam.

v) Recording system →

The signal from the detector is finally received by the recording system, it can be done by a recorder pen.

UV - Visible Spectrophotometer

1) Single beam UV spectrophotometer :-



• Step involved :-

i) UV radiation is given off by the

source

ii) A convex lens gathers the beam of radiation

and focuses it on the slit.

iii) The inlet slit permits light from the source to pass, but block-out stray radiation

iv) The light then reaches the monochromator, which split it up according to wavelength.

v) The exit slit is positioned to allow light of the required wavelength to pass through.


vi) Radiation at all other wavelength is blocked out

vii) The selected radiation passes through the sample cell to the detector, which measures the intensity of the radiation reaching it.

viii) By comparing the intensity of radiation before and after it passes through the sample, it is possible to measure how much radiation is absorbed by the sample at the particular wavelength used.

ix) The output of the detector is usually recorded on graph paper.

Double Beam Spectrophotometer



- A tungsten lamp (400-800 nm)
- A D₂ lamp (200-400 nm)

• Step involved:—

i) The radiation from the selected source passes through a fixed slit to the surface of reflection grating (monochromator), and from the diffracted radiation, the desired wavelength is selected.

ii) The selected beam of light falls on V-shaped mirror, called a beam splitter that splits the radiation into two beams; one of which passes through the reference cell containing pure solvent and other simultaneously passes through the sample solution cell.

iii) The transmitted light from the two cells go to the photoelectric detector alternatively, and the difference in the absorbance by the solvent and the sample solution is measured electronically with very high accuracy.

• Applications of UV Visible Spectroscopy

UV-Vis spectroscopy has been mainly applied for :-



i) Detection of functional groups - It is applied to detect the presence

or absence of a chromophore.

ii) Extent of conjugation - conjugation can

shifts the absorption to longer wavelengths. (\uparrow absorption \neq \uparrow wavelength)

iii) Identification of a unknown ~~(a)~~ ~~(b)~~

Compound - By comparing the spectrum of an unknown compound with that of a known compound.

iv) Detection of impurities - Additional peak can be due to impurities in sample

and can be compared with that of standard materials.

v) Difference in conjugated & non-conjugated compounds - By comparing their wavelength.

* Quantitative analysis :-

i) Spectrophotometric titrations

ii) Single component analysis

iii) Multi component analysis

i) Spectrophotometric titration \rightarrow



It is the process of



determining the quantity of a sample by adding measured increments (increased)

of a titrant until the end point...

The endpoint is where the graph is discontinuous.

- The titrations are based on Beer's law $A \propto C$
- Titration curve is plot of absorbance Vs volume of titrant

Titration curves



ii) Single component analysis →

- If sample is consist of only drug that absorbs radiation and excipient does not absorb radiation at λ_{max} of drug. λ_{max} is the wavelength at which a substance absorbs maximum radiation.
- Methods use - Standard absorbability value, By using $A = abc$, By using Beer's curve. etc.

iii) Multicomponent analysis →

$$A \propto C$$

If some sample has more than one drug and all of them absorb at same wavelength then this analysis is required

[or]

If sample carry excipients, impurities product which absorb at λ_{max} of drug then multicomponent analysis is required.

• Methods used -

- Simultaneous equation method
 - Difference spectrophotometry.
 - Absorptivity factor method
 - Dual wavelength method
- etc...

XX — Completed — XX

MUST WATCH VIDEOS ON

"CAREWELL PHARMA" YOUTUBE

CHANNEL for BETTER
UNDERSTANDING.

Follow us on :-

YOUTUBE

TELEGRAM

INSTAGRAM